

## 1. Unpacking and Setup

- **Unbox** the Rhythm P2 Portable Oxygen Concentrator carefully.
- **Check contents:** Ensure the device, AC and DC power cords, battery, carrying case, nasal cannula, and user manual are included.
- **Charge the battery:** Plug the AC or DC adapter into the device and charge the battery fully before first use.

## 2. Powering the Device

- **Attach the battery:** Align and snap the battery into place at the bottom or back of the unit.
- **Power on:** Press the power button (usually at the top/front) to turn the device on.
- **Battery status:** Check the battery level on the display screen to ensure enough charge for use.

## 3. Connecting the Cannula

- **Attach the nasal cannula:** Insert one end of the nasal cannula tubing into the oxygen outlet port on the side or front of the device.
- **Place the cannula:** Put the prongs of the nasal cannula into your nostrils and loop the tubing over your ears to secure it.

## 4. Adjusting Oxygen Flow

- **Set the flow rate:** Use the buttons or knob to adjust the oxygen flow to your prescribed setting. Flow rate is shown on the display screen.
- **Monitor output:** Ensure oxygen is flowing properly by feeling the airflow through the nasal cannula.

## 5. Using the Device

- **Portable use:** Once powered on, you can carry the device using the provided strap or handle.
- **Battery monitoring:** Regularly check the battery status on the screen to ensure you don't run out of power during use. Charge when necessary.

## 6. Charging the Battery

- **Charge at home:** Plug the AC adapter into a wall socket and the other end into the concentrator to charge.
- **Charge in the car:** Use the DC power cord to charge the concentrator through your vehicle's power outlet.

## 7. Maintenance

- **Clean the exterior:** Wipe down the outside of the device with a soft, dry cloth.
- **Cannula care:** Replace the nasal cannula regularly as recommended by your healthcare provider.
- **Check filters:** Periodically inspect the air intake filter and replace it as specified in the manual.

## 8. Powering Off

- **Turn off the device:** Press and hold the power button until the device shuts down.
- **Store safely:** When not in use, store the P2 in a cool, dry place.

## Troubleshooting

- **Low battery warning:** Charge the battery if you hear a low battery alarm or see a low battery warning on the display.
- **No oxygen flow:** Ensure the cannula is properly connected, the device is powered on, and the flow setting is correctly adjusted.

## Safety Tips

- Do **not** block the air intake vent.
- Use **only** approved power cords and accessories.
- Avoid **smoking** or using the device near open flames.

For any technical issues or detailed maintenance, refer to the user manual or contact the manufacturer.